COVID-19 (SARS-CoV-2) Antigen Test Kit (Colloidal Gold) Saliva


For self-testing. For in vitro diagnostic use only. Please read the instruction carefully before use.


Intended use

This product is used for in vitro qualitative detection of the SARS-CoV-2 antigen in human saliva specimen. This product is intended for home self-testing as a rapid test for novel coronavirus infection. Both symptomatic and asymptomatic infections can be tested. However, please do not make a medical relevance decision without consulting with your doctor. This product is suitable for users aged 10 and over. Parents or eligible adults should assist users under the age of 10 to complete the test.

Materials and Components

Preparation before the test

Interpretation of test results

Negative result: If there is only a control line (C) and the test line (T) is colorless, it indicates that SARS-CoV-2 antigen has not been detected and the result is negative. If the test result is positive, it indicates that SARS-CoV-2 antigen has been detected and the result is positive. **Currently, there is a suspected infection of COVID-19. Contact your doctor or local health department immediately. Comply with local regulations, self-isolate and report according to local regulations. Perform PCR test for confirmation.**

Quality Control

Program control is included in the test. A red line appearing in the control region (C) is the internal procedural control. It confirms sufficient volume of the specimen.

Performance Index

- Limit of detection (LOD):** T_{CD₅₀} is 80.
- High Dose Hook Effect:** When the virus concentration exceeds 1.4 x 10⁷ TCID₅₀/mL, the result may be false negative.
- Cross-reactivity:** There is no cross-reactivity, including human coronavirus 229E, human coronavirus OC43, human coronavirus NL63, human coronavirus HKU1, MERS-coronavirus, SARS coronavirus, adenovirus 3, and parainfluenza virus type 2, Enterovirus, respiratory syncytial virus (A), parainfluenza virus type 3, parainfluenza virus type 4a, Influenza A H3N2 (Wisconsin/67/05), Influenza A H1N1, Influenza B (VICTORIA), Rhinovirus (HRV430), Haemophilus influenzae Streptococcus pneumoniae, Streptococcus pyogenes, Candida albicans, Bacillus pertussis, Mycoplasma pneumoniae, Legionella pneumophila, Staphylococcus aureus, Streptococcus pneumoniae, Mycoplasma pneumoniae, Chlamydia pneumoniae, Legionella pneumophila, Mycobacterium tuberculosis, Pneumocystis, Pseudomonas, Bacteria human pneumonia virus (hMPV), parainfluenza virus type 1, Staphylococcus epidermidis, Streptococcus salivarius, etc.
- Microbial Interference Studies:** There is no interference in studies on the following microorganisms, including parainfluenza virus type 1, parainfluenza virus type 2, parainfluenza virus type 3, parainfluenza virus type 4a, adenovirus, human pneumonia virus (hMPV), A H3N2 Influenza (Wisconsin/67/05), H1N1 influenza, Haemophilus influenzae, Streptococcus pneumoniae, Streptococcus pyogenes, Influenza B (Malaysia/2506/04), enterovirus, respiratory syncytial virus, Rhinovirus, Chlamydia pneumoniae, Legionella pneumophila, Mycobacterium tuberculosis, Pneumocystis, Pseudomonas, Candida albicans, Bacillus pertussis, Mycoplasma pneumoniae, Staphylococcus epidermidis, Streptococcus salivarius, human coronavirus 229E, human coronavirus OC43, human coronavirus NL63, human coronavirus HKU1, MERS coronavirus, etc.
- Interference Studies:** There is no interference in studies on the following substances, including blood, mucin, Alkaloid, dexamethasone, Nelmed, benzocaine, tobacramin, mucipronin, biotin, etc.

Clinical Performance

The overall study scale was 800 cases, 150 positive samples and 450 negative samples. Statistics of test results of saliva samples:

Reference RT-PCR Assay		95% Wilson Score CI	
	LCI	UCI	
DEEPBLUE SARS-CoV-2 ANTIBODY TEST KIT	POS 148	0	148
	NEG 2	450	452
	TOTAL	150	450
	PPA	>99.9%	98.17%
	NPV	>99.9%	98.17%
	PPV	>99.9%	98.17%
	NPV	>99.9%	98.17%

Sensitivity: 98.7% (95% CI: 92.33% - 99.07%)
Specificity: >99.9% (95% CI: 98.17% - 100%)
Sensitivity: Compared with the RT-PCR Assay, among people infected with SARS-CoV-2 virus, the probability of correct detection by the DeepBlue SARS-CoV-2 Ag Test Kit.
Specificity: Compared with the RT-PCR Assay, among people who have not been infected with SARS-CoV-2 virus, the probability of correct detection by the DeepBlue SARS-CoV-2 Ag Test Kit.

Test principle

This product uses the double antibody sandwich method to detect the SARS-CoV-2 N protein. When the sample contains the coronavirus antigen, both the test line (T) and the control line (C) will appear, and the result will be positive. When the sample does not contain the coronavirus antigen or no coronavirus antigen is detected, the test line (T) will not appear, only control line (C) will appear.

Limitations of inspection methods

- This test kit is only used for in vitro diagnosis.
- This test kit is only used to detect human saliva. The results of other specimens may be wrong.
- This test kit is only used for qualitative detection and cannot indicate the level of SARS-CoV-2 antigen in the specimen.
- This test kit is only a clinical auxiliary diagnostic tool. If the result is positive, it is recommended to use other methods for further examination in time and the doctor's diagnosis shall prevail.
- This test does not determine the etiology of the respiratory infection caused by microorganisms other than the SARS-CoV-2 virus.
- This test can detect both the viable and the non-viable SARS-CoV-2 virus. The accuracy of the test depends on the quality of the swab sample—false negative results may be given following poor sampling.
- Any failure to respect the test procedure may negatively impact the performance of the test and/or invalidate the test result.
- If the result of the test is negative, vital symptoms persist, it is advised that you carry out additional tests using other clinical methods. A negative result at no time rules out the presence of antigens of the SARS-CoV-2 virus in the sample, as they may be present but at a level inferior to the minimum detection level of the test, or if the sample has been collected incorrectly.
- A negative result does not rule out infection by the SARS-CoV-2 virus, particularly in people who have come into contact with the virus. Follow-up tests with molecular diagnostics should be scheduled to rule out infection in these people. Persons who show symptoms of the disease but have a negative result until infection is ruled out should follow current restrictions.
- This test is not a substitute for a medical consultation, or for the result of a biological analysis carried out in a medical analysis laboratory.
- Positive test results do not exclude the possibility of co-infections of other pathogens.

Warnings and Precautions

- Read the instructions carefully before using the kit and strictly control the reaction time. If you do not follow the instructions, you will get inaccurate results.
- Do not eat, drink, chew gum, smoke or vape for at least 30 minutes before collecting saliva. False negative results can occur if the saliva is not collected properly.

- Guard against moisture, do not open the aluminum platinum bag before it is ready for testing. Do not use it if the aluminum foil bag is damaged or the test device is damp. Please use it within the validity period.
- Do not replace the components in this kit with components in other kits.
- Do not dilute the specimen when tested, otherwise you may get inaccurate results. The kit shall be stored in strict accordance with the conditions specified in this manual. Please do not store the kit under freezing conditions.
- The test methods and results must be interpreted in strict accordance with this specification.
- Negative results may occur if the SARS-CoV-2 antigen titer in the specimen falls below the minimum detection limit of this kit.
- If the extraction reagent is individual packaging and one piece per test device, the batch number, expiration date and other information cannot be marked separately due to the space is limited, but those information will be consistent with the corresponding test kit.
- There is no reduction in sensitivity in the DeepBlue Antigen test against the UK variant of SARS-CoV-2.
- Do not swallow the extraction reagent. If accidentally touch the human skin, eyes or mucous membranes, please rinse with water immediately. If discomfort occurs, please consult a doctor.

Storage conditions & period of validity

- Store at 4-30°C, and it is valid for 36 months.
- After the aluminum foil bag is unsealed, the test device should be used as soon as possible and within one hour (15 ~ 30°C, Humidity ≤ 80%).

Sample Transport and Storage

Freshly collected specimens should be processed as soon as possible. It should be no later than one hour after collection.

Quality Control

Program control is included in the test. A red line appearing in the control region (C) is the internal procedural control. It confirms sufficient volume of the specimen.

Performance Index

- Limit of detection (LOD):** T_{CD₅₀} is 80.
- High Dose Hook Effect:** When the virus concentration exceeds 1.4 x 10⁷ TCID₅₀/mL, the result may be false negative.
- Cross-reactivity:** There is no cross-reactivity, including human coronavirus 229E, human coronavirus OC43, human coronavirus NL63, human coronavirus HKU1, MERS-coronavirus, SARS coronavirus, adenovirus 3, and parainfluenza virus type 2, Enterovirus, respiratory syncytial virus (A), parainfluenza virus type 3, parainfluenza virus type 4a, Influenza A H3N2 (Wisconsin/67/05), Influenza A H1N1, Influenza B (VICTORIA), Rhinovirus (HRV430), Haemophilus influenzae Streptococcus pneumoniae, Streptococcus pyogenes, Candida albicans, Bacillus pertussis, Mycoplasma pneumoniae, Legionella pneumophila, Staphylococcus aureus, Streptococcus pneumoniae, Mycoplasma pneumoniae, Chlamydia pneumoniae, Legionella pneumophila, Mycobacterium tuberculosis, Pneumocystis, Pseudomonas, Bacteria human pneumonia virus (hMPV), parainfluenza virus type 1, Staphylococcus epidermidis, Streptococcus salivarius, etc.
- Microbial Interference Studies:** There is no interference in studies on the following microorganisms, including parainfluenza virus type 1, parainfluenza virus type 2, parainfluenza virus type 3, parainfluenza virus type 4a, adenovirus, human pneumonia virus (hMPV), A H3N2 Influenza (Wisconsin/67/05), H1N1 influenza, Haemophilus influenzae, Streptococcus pneumoniae, Streptococcus pyogenes, Influenza B (Malaysia/2506/04), enterovirus, respiratory syncytial virus, Rhinovirus, Chlamydia pneumoniae, Legionella pneumophila, Mycobacterium tuberculosis, Pneumocystis, Pseudomonas, Candida albicans, Bacillus pertussis, Mycoplasma pneumoniae, Staphylococcus epidermidis, Streptococcus salivarius, human coronavirus 229E, human coronavirus OC43, human coronavirus NL63, human coronavirus HKU1, MERS coronavirus, etc.
- Interference Studies:** There is no interference in studies on the following substances, including blood, mucin, Alkaloid, dexamethasone, Nelmed, benzocaine, tobacramin, mucipronin, biotin, etc.

Clinical Performance

The overall study scale was 800 cases, 150 positive samples and 450 negative samples. Statistics of test results of saliva samples:

Reference RT-PCR Assay		95% Wilson Score CI	
	LCI	UCI	
DEEPBLUE SARS-CoV-2 ANTIBODY TEST KIT	POS 148	0	148
	NEG 2	450	452
	TOTAL	150	450
	PPA	>99.9%	98.17%
	NPV	>99.9%	98.17%
	PPV	>99.9%	98.17%
	NPV	>99.9%	98.17%

Sensitivity: 98.7% (95% CI: 92.33% - 99.07%)
Specificity: >99.9% (95% CI: 98.17% - 100%)
Sensitivity: Compared with the RT-PCR Assay, among people infected with SARS-CoV-2 virus, the probability of correct detection by the DeepBlue SARS-CoV-2 Ag Test Kit.
Specificity: Compared with the RT-PCR Assay, among people who have not been infected with SARS-CoV-2 virus, the probability of correct detection by the DeepBlue SARS-CoV-2 Ag Test Kit.

Index of Symbols

	In vitro diagnostic medical device
	Use-by date
	Temperature limit
	Authorized representative in the European Union
	Do not reuse
	Consult instructions for use
	Manufacturer
	Unique device identifier

	Batch code
	Catalogue number
	Keep dry
	Keep away from sunlight
	Do not use if package is damaged and consult instructions for use
	Contains sufficient for n tests
	CE Mark
	Importer
	Sterilized swab

Shenzhen KangDaAn Biological Technology Co., LTD. East-1, 3rd floor, Building 2, Shunheida factory, Lixiangdong industrial zone, Xili street, Nanshan district, Shenzhen China

Share Info Consultant Service LLC Repräsentanzbüro Herdorf, Lohweg 63, 40549 Düsseldorf, Germany

Goodwood Medical Care Ltd. 1-2 Floor, 3-919 Yongzheng Street, Jinzhou District, Dalian, China

CMC Medical Devices & Drugs S.L. C/ Horacio Lengua No° 18, CP 29006, Málaga, Spain

Hangzhou Yigouren Biotechnology Co., Ltd. Room 201, Building 5, No.360, Tianshuang Road, Yuhang Street, Yuhang District, Hangzhou, 311100 Zhejiang, P.R.China.

BioDomi GmbH Schokkalandring, 13 65795 Hattersheim am Main, Germany

ANHUI DEEPBLUE MEDICAL TECHNOLOGY CO., LTD. No. 777 Jingshan Road, High-Tech Development Zone 230088 Hefei, Anhui, PEOPLE'S REPUBLIC OF CHINA

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Shenzhen KangDaAn Biological Technology co., LTD, East-1, 3rd floor, Building 2, Shungheda factory, Luxiangdong industrial zone, Xili street, Nanshan district, Shenzhen China

Share Info Consultant Service LLC Repräsentanzbüro Heerdter Lohweg 83, 40549 Düsseldorf, Germany

Goodwood Medical Care Ltd, 1-2 Floor, 3-619 Yongzheng Street Jinhzhou District, Dalian, China

CMC Medical Devices & Drugs S.L, C/ Horacio Lengua No° 18, CP 29006, Málaga, Spain

Hangzhou Yiguoen Biotechnology Co., Ltd, Room 201, Building 5, No.360, Tianmushan Road, Yuhang Street, Yuhang District, Hangzhou, 311100 Zhejiang, P.R.China,

BioDom GmbH Schokoladenring,13 65795 Hattersheim am Main, Germany

COVID-19 (SARS-CoV-2) Antigen Test Kit (Colloidal Gold) - Saliva

ANHUI DEEPBLUE MEDICAL TECHNOLOGY CO., LTD, No. 777 Jimgingshan Road, High-Tech Development Zone 230088 Hefei, Anhui, PEOPLE'S REPUBLIC OF CHINA

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Luxus Lebenswelt GmbH Kochstr. 1, 47877, Willich, Germany