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
ACCESSORI

Accessorio	Produttore	Representante CE	Marcio CE
Tampone A	Goodwood Medical Care Ltd.	CMC Medical Devices & Drugs S.L. C/Horacio Lengo No.18,CP29006,Málaga,Spain	CE 0197 acc.93/42/CEE
Tampone B	Jiangsu Rongye Technology Co.,Ltd Toujiao Town, Yangzhou City, 225109 Jiangsu, P.R.China	Riomavix S.L. Calle de Almansa 55,1D, Madrid 28039 Spain	CE 0197 acc.93/42/CEE
Tampone C	Dalian Rongbang Medical Healthy Devices Co.,Ltd. Macyingzi Hamlet, Dalianwan Town, Ganjingzi District, Dalian, 116113	Lotus NL B.V. Koningin Julianaplein10, le Verd,2595AA,The Hague,Netherlands	CE 0197 acc.93/42/CEE

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INTENDED USE

The COVID-19 Antigen Rapid Test Kit (Swab) is a single-use test kit intended to detect the novel coronavirus SARS-CoV-2 that causes COVID-19. This test is intended for home use with self-collected nasal swab samples in individuals aged 18 and older. Sampling from anyone under the age of 18 should be performed under the guidance of an adult. People who are unable to carry out the test on their own should seek support. This test utilizes the lateral flow method of nucleocapsid protein antigen in individuals with suspected COVID-19. Positive results are indicative of the presence of SARS-CoV-2. Individuals who test positive should self-isolate and seek additional care from their healthcare provider. Positive results do not rule out bacterial infection or coinfection with other viruses. Negative results do not preclude SARS-CoV-2 infection. Individuals who test negative and continue to experience COVID-like symptoms should seek follow up care from their healthcare provider.

PRINCIPLE

The COVID-19 Antigen Rapid Test Kit (Swab) is an immunochromatographic membrane assay that uses highly sensitive antibodies to detect SARS-CoV-2 nucleocapsid protein from nasal swab specimens. Anti-SARS-CoV-2 specific antibodies are immobilized onto the test region of the membrane and combined with other reagents/pads to construct a test strip. During testing, the specimen reacts with anti-SARS-CoV-2 antibodies conjugated to colored particles and precolated onto the sample pad of the test. The mixture then migrates upward on the membrane chromatographically by capillary action and reacts with the reagents in test line region. Therefore, if the specimen contains N antigen of SARS-CoV-2, a colored line will appear in test line. If the specimen does not contain N antigen of SARS-CoV-2, no colored line will appear in the test line regions, indicating a negative result. To serve as a procedural control, a colored line will always appear in the control line region, indicating that the proper volume of specimen has been added and membrane wicking has occurred. The test is designed to detectation of nucleocapsid protein antigen in nasal swab specimens.

KIT COMPONENTS

Component	1 Test /Kit	5 Tests /Kit	10 Tests /Kit	20 Tests /Kit	25 Tests /Kit
COVID-19 Antigen Test	1	5	10	20	25
Extraction tube with buffer	1	5	10	20	25
Sterilized nasal swab	1	5	10	20	25
Waste bag	1	5	10	20	25
Workstation	1	1	1	1	1
Package insert	1	1	1	1	1

ADDITIONAL SPECIAL EQUIPMENT

Timer

WARNINGS AND PRECAUTIONS

- Do not use after expiration date. Do not use if kit is damaged or open. Do not reuse the tests. Do not mix components from different kit lots.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Handle all specimens as if they contain infectious agents. Discard the using testing materials in accordance with local regulations.
- The extraction buffer contains a salt solution that the solution contacts the skin or eyes. Flush with copious amounts of water. Don't swallow the buffer. When swallowing the buffer, rinse the mouth thoroughly with water and give plenty of water to dilute the substance. If any discomfort, seek medical attention immediately.
- Children and elder please use the test under the guardian.

STORAGE AND STABILITY

Store unused test devices unopened at 4°C-30°C. Do not freeze the kit or expose the kit over 30°C. If stored at 4°C-6°C, ensure that the test device is brought to room temperature before opening. The shelf life of the test device is 36 months. The detailed expiration date printed on the sealed pouch and outer box.

TEST PROCEDURE

Open the kit box. Check the components before use.

Please read all instructions carefully before you begin.

[Preparation before sampling]

- Get a flat area ready, like a table. Make sure it is clear, clean and dry.
- Take off hand jewellery.
- Wash your hands for 20 seconds. Use soap and warm water, or hand sanitizer. Dry your hands using clean disposable paper towels.


For better protection and avoid cross-contamination, disposable gloves, masks and eye protection (not provided in the package) are recommended.

[Specimen Preparation]

- Press the hole on the box, peel off the foil film on the extraction tube and insert into the hole. (For 1 test/kit and 5 tests/kit, insert the extraction tube into the hole of the box).
- Open nasal swab package at the sticky end and take the nasal swab out.
- Insert the soft tip into left nostril and insert into the nostril 2.5cm (1 inch) for adult. Note: For child, the maximum depth of insertion into the nostril may be less than 2.5 cm and should be carefully and appropriately adjusted by the person, who collects sample.
- Firmly brush against the inside of the nostril in a circular motion 5 times or more.
- Move the nasal swab to the right nostril and repeat the previous action. Make sure an adequate sample is collected.
- Insert the nasal swab into the tube which contains extraction buffer.
- Rotate nasal swab at least 8-10 times while pressing nasal swab against the bottom of the side of the tube of the tube.
- Insert the nasal swab into left nostril and insert into the tube to release as much liquid as possible.
- Cover the tube with dropper tip tightly and insert the tube back into the hole.

[Test Procedure]

- Open the sealed pouch and take out the test cassette. For best results, the test should be performed in one hour.
- Hold the tube vertically upside down over the sample well.
- Add 3 drops specimen into the sample well by gently squeezing the sides of the tube, then start the timer

4. Wait for colored lines to appear. The test result can be read in 10-15 minutes. DO NOT read after 20 minutes.

[After the testing]

- After you have done the test, put all parts of the kit in the waste bag. Discard the waste bag in accordance with local regulations.
- Clean the table and timer button with alcohol or surface sanitizer. Wash your hands after each test.

INTERPRETATION OF TEST RESULT

POSITIVE RESULT

 A colored line appears in the control line region (C) and a colored line appears in test line region (T).
***NOTE:** The intensity of the color in the test line region will vary depending on the concentration of N antigen of SARS-CoV-2 in the specimen. Therefore, any shade of color in the test line region should be considered positive.

NEGATIVE RESULT

 A colored line appears in the control line region (C) and no line appears in test line region (T).

INVALID RESULT

 No line appears in control line region (C). Insufficient buffer volume or incorrect procedural techniques are the most likely reasons for control-line failure. Review the procedure and repeat the procedure with a new test device. If the problem persists, discontinue using the test kit immediately and contact your distributor.

BUILT-IN CONTROL

This test contains a built-in control feature, the C line on the test. The C line develops after adding sample solution. Otherwise, review the whole procedure and repeat test with a new device.

WHAT SHOULD I DO AFTER TEST

If the test result is positive	A positive test result means that SARS-CoV-2 virus was detected in your sample and it is very likely you have COVID-19. <ul style="list-style-type: none"> Contact your doctor / general practitioner or the local health department immediately Comply with local guidelines on self-isolation.
If the test result is negative	A negative result means that the SARS-CoV-2 virus was not detected in the sample, so COVID-19 is not confirmed; however, the test cannot exclude this option completely, and additional testing may be needed. <ul style="list-style-type: none"> If it is suspected, repeat the test after 1 - 2 days, as the coronavirus cannot be accurately detected in all phases of an infection. If you test negative and continue to experience COVID-19-like symptoms, you should seek follow up care with your health care provider.
If the test result is invalid	Possibly caused by incorrect test execution <ul style="list-style-type: none"> Repeat the test If the test results remain invalid, contact a doctor or a COVID-19 test center.

Note: Do not take any decision of medical relevance without first consulting your medical practitioner.

PERFORMANCE CHARACTERISTICS

1. Limit of Detection: The Limit of Detection (LoD) was determined by evaluating different dilutions of heat inactivated SARS-CoV-2 virus. LoD is defined as the virus concentration 1.3×10^4 TCID₅₀/mL at which a minimum of 19 replicates out of 20 generate a Reactive result.
 2. High-dose Hook Effect: No hook-effect was observed at 3.33×10^7 TCID₅₀/mL of inactivated SARS-CoV-2 virus culture.
 3. Clinical study: A comparison was conducted with RT-PCR, according the test data listed below table. The Relative Sensitivity is 96.4% (Total 220 RT-PCR (nasopharyngeal swab) confirmed positive samples: 212 RT-PCR (nasopharyngeal swab) confirmed positive samples were correctly detected by the COVID-19 rapid antigen test. There are 8 false negative cases.), The Relative Specificity is >99.9 (Total 200RT-PCR(nasopharyngeal swab)-confirmed negative samples: 200 RT-PCR(nasopharyngeal swab)-confirmed negative samples were correctly detected by COVID-19 antigen rapid test. There are no false positive) and the overall agreement is 98.1% (Total 420RT-PCR (nasopharyngeal swab) confirmed samples: 412 RT-PCR (nasopharyngeal swab) confirmed samples were correctly detected by the COVID-19 rapid antigen test).

		PCR result		Total
		Positive	Negative	
Acleu Medical Test	Positive	212	0	212
	Negative	8	200	208
	Total	220	200	420

4. Cross-reactivity: Cross-reactivity studies are performed to demonstrate that the test does not react with the following microorganisms in the table below at concentration of 1×10^4 TCID₅₀/mL for viruses and 1×10^4 CFU/ml for bacteria.

Human metapneumovirus (hMPV)	Human parainfluenza virus 1	Rhinovirus	<i>Bordetella pertussis</i>	<i>Streptococcus pneumoniae</i>
Human coronavirus OC43	Human parainfluenza virus 2	Enterovirus	<i>Chlamydia pneumoniae</i>	<i>Streptococcus pyogenes</i>
Human coronavirus 229E	Human parainfluenza virus 3	Adenovirus	<i>Haemophilus influenzae</i>	<i>Mycobacterium tuberculosis</i>
Human coronavirus NL63	Human parainfluenza virus 4	Influenza A	<i>Legionella pneumophila</i>	<i>Staphylococcus aureus</i>
Respiratory Syncytial Virus	MERS	Influenza B	<i>Mycoplasma pneumoniae</i>	<i>Candida albicans</i>
Human coronavirus HKU1				

5. Interference: The following endogenous interference substances were evaluated at the concentrations listed and no effect was found. Whole blood (2%), three OTC nasal sprays (10%), three OTC nasal drop (25%), three nasal mouthwashes (25%), 4-Acetamidophenol (10mg/ml), Acetylsalicylic acid (20mg/ml), Chlorpheniramine (5 mg/ml), Dextromethorphan (10mg/ml), Diphenhydramin (5mg/ml), Ephedrine (20mg/ml), Guaiacol glyceryl ether (20mg/ml), Oxymetazoline (10mg/ml), Phenylephrine (100mg/ml), Phenylpropanolamine (20mg/ml), Oseltamivir Phosphate (10mg/ml), Mupirocin (10mg/ml), Vitamin A (10%), D-Panthenol (10%)

LIMITATIONS AND POSSIBLE ERRORS

- The COVID-19 Antigen Rapid Test Kit (Swab) is intended for use as a self-test and may only be used for qualitative detection of SARS-CoV-2 antigen. The colour intensity of a positive line shall not be evaluated as quantitative or semi-quantitative.
- The COVID-19 Antigen Rapid Test Kit (Swab) should be only for the detection of SARS-CoV-2 antigen, not for any other viruses or pathogens.
- The performance was evaluated using the procedures provided in this product insert only. Modifications to these procedures may alter the performance of the test.
- A negative result does not completely exclude the possibility of SARS-CoV-2 infection.
- The results obtained with this assay, especially in the case of weak test lines which are difficult to interpret, should retest or go to a medical institution for testing.
- The test is intended for infection detection and not for determine infection status. The test is used for the auxiliary diagnosis of patients with the COVID-19 and can not be used as the sole diagnostic indicator of whether the test subject is infected with the COVID-19.

ACCESSORIES

Accessory	Manufacturer	EC Representative	CE Mark
Swab A	Goodwood Medical Care Ltd, 1-2 Floor,3-919 Yonzheng Street,Jinzhou District Dalian, Liaoning China,116100	CMC Medical Devices & Drugs S.L./CHoracio Lengo No.18,CP29006,Málaga,Spain	
Swab B	Jiangsu Rongye Technology Co.,Ltd Toudiao Town,Yangzhou City,225109 Jiangsu, P.R.China	Riomavix S.L. Calle de Almansa 55,1D, Madrid 28039 Spain	
Swab C	Dalian Rongbang Medical Healthy Devices Co.,Ltd, Maoyingzi Hamlet, Dalianwan Town, Ganjingzi District, Dalian, 116113	Lotus NL B.V. Konigin Julianaplein10, le Verd,2595AA, The Hague,Netherlands	

INDEX OF SYMBOLS

	Do not reuse		For in vitro diagnostic use only
	Stored between 4-30 °C		Consult instruction for use
	Keep away from sunlight		Lot number
	Use-by date		Contains sufficient for <-> tests
	Manufacturer		Date of Manufacture
	Sterilized using ethylene oxide		Sterilized using irradiation
	Catalogue Number		CE Mark
	Authorized Representative in the European Community		NIC GmbH Erlenweg 13, 49076 Osnabrück, Germany Tel.: +49 541 911 67 06 E-mail: info@nic-industry.com

POSITIVE RESULT

 A colored line appears in the control line region (C) and a colored line appears in test line region (T).
***NOTE:** The intensity of the color in the test line region will vary depending on the concentration of N antigen of SARS-CoV-2 in the specimen. Therefore, any shade of color in the test line region should be considered positive.

NEGATIVE RESULT

 A colored line appears in the control line region (C) and no line appears in test line region (T).

INVALID RESULT

 No line appears in control line region (C). Insufficient buffer volume or incorrect procedural techniques are the most likely reasons for control-line failure. Review the procedure and repeat the procedure with a new test device. If the problem persists, discontinue using the test kit immediately and contact your distributor.

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